

## How to restart the cooling of the crystallization cold room i.e. after a power cut.

Enter the crystallization room (VF 00.015) and get the key of the technic cabinet VF 00.020. It is hanging, at position 20, in the key box situated on the wall on the left hand side when you enter the crystallization room.



Use the key to open the door of the technic cabinet VF 00.020. The cabinet is one on your left hand side, when you leave the crystallization room. Next to the room number, there is a “Kälte” label. Open the door with the key. You should see the following panel:



In case of problem, the red light should be on saying “Störung”. If that’s the case, and it usually is, after power cut, press the green button to restart the machine. It usually makes some noise when running, so you should hear it when it restarts. This allows circulation of water cooled on the rooftop.

Close, and lock, the door and hang it back it the key box of the crystallization room.

Then go “behind” the crystallization room, along the windows. There you’ll find the cooling device that uses the cold water circulating in the pipe to cool the coldroom.



In front of it, on the small metal box (lower right of the image), you’ll find a “reset” button. Press it once to restart that device as well.

In principle, you should see that the temperature of the cold room will decrease quickly after the restart of both devices